


SUPPLEMENTARY INFORMATION

Schematics, bill of materials, notes on PCB impression and on the first use.

SOMA







Figure 2. View of PCB card (components side).

List of SOMA components (with LCSC.com references) :

Cm1 = 10 μ F (Cm1 \geq CF)	C162684	C531341	50Vdc 25Vdc
CF = 10 μ F	C162684	C531341	
CS = 4,7 μ F	C2839238	C3696971	50Vdc 25Vdc
RSB = 4k7		C144486	
RS = 5K6		C229724	
RF = — (Not used)			
Rin1 = 10K		C140407	
Rm = 1M		C108083	
Rg = 1K		C102060	
RA = RB = 100K		C144522	
RS = 5K6		C229724	
RP1 = 4k7		C144486	
RP2 = 1k2		C212484	
RP4 = — (not used)			
P1 = 2K (multiturn potentiometer)		C116322	
R7 = 0 (short-circuit)			
R10 = 10K		C140407	
D3 is replaced with résistance RD = 1K5			
R14//R15 = (47k//100k) \pm 33K (to be adjusted)	C144490	// C144522	
RP3//RP6 = \pm 1K5 (adjustment of the peak voltage with <u>TEA</u> to +4Vdc)		C114929	
RREP1//RREP2 = 1K5//2K7 (adjustment of the plateau voltage with <u>TEA</u> to +2Vdc)	C114929	//	
C104725			
D1 = D2 = DF = 1N4148			
R2 = R3 = 180K		C104659	
R4 = 47K		C102201	
R5 = 10K		C140407	
R6 = 47K		C102201	
R13 = 47K		C102201	
RA1 = 12K		C171051	
RAL1 et RAL2 = 1K		C173143	
R1 = 47K		C102201	
DF, D1, D2 = 1N4148		C402212	
LDR = (100K à 200K), P=3mm		C125631	
U1 et U2 = LM321TR		C2842352	
TR1 à TR3 = 2N3904		C8667	
DZ1 = 12Vdc		C388130	
Contact pins = RH-5011		C5199814	
Male connector HDX1 3P (Power 24Vdc)	C441172	KF2EDGR-3.5-3P	
Female connector HDX2 8P (DENDRITE)	C2897390		
Female connector HDX3 5P (Axon Hill)	C35167		
Female connector HDX4 10P (extension cards)	C2932683		
1 WHITE LED	C965818		
1 BLUE LED	C84259		
3 RED LEDs	C965812		
SWITCH double-inverter =	C189615		

NODE OF RANVIER



List of Nove of Ranvier components (with LCSC.com references) :

CX1=4,7 μ F // CX13 =2,2 μ F = 6,9 μ F ceramic	C162274 // C237392
CX2=4,7 μ F	C162274
CX6 = 47nF	C107154
RX4=5,6K	C229724
RX9 = 1M	C108083
RX8 = 1K2	C212484
RX6 = 4K7	C144486
RX26 = 100	C245445
RX28 = 4K7	C144486
RX2 = RX3 = 100K	C144522
RX1 = 1K + PX2 (regulation of signal delay)	C102060 + C116322
RXSHUNT = 1M	C108083
RX60 = 33K	C137309
RX61 = 2k2	C114933
RX32 = 56k	C137251
RX18 = 12K	C171051
RX31 = 33K	C137309
RXgain = 330	C17930
PX1 and PX2 = 2K (multi-turn potentiometers)	C116322
TR1 et TR5 = 2N3904	C8667
UX1 à UX3 = LM321tr	C2842352
1 WHITE LED (Signal)	C965818
1 BLUE LED (Ionic Gradients)	C84259
Contact pins	C5199814
Female connector 5P (output)	C35167
Male connector 5P (input)	C138801

The following components on the PCB card are not used

CX12 = ---
CX15 = ---
CX14 = ---
CX10 = ---
CX5 = ---
CX3 = — (in // to CX2)

DENDRITES




Figure 4. Schematics of dendrites




Figure 5. View of PCB card of dendrites (components side).

List of dendrites components (with LCSC.com references):

PX1 = 2k (multi-turn potentiometer)	C116322
RX30 = 1K2	C17891
RX27 = —	
RX1, RX8 = 82K	C17979 C17979
RX2, RX9 = 82K	C17979 C17979
RX3, RX10 = 82K	C17979 C17979
RX5, RX12 = 150K	C2907431 C114931
RX6, RX13 = 150K	C2907431 C114931
RX7, RX14 = 150K	C2907431 C114931
RX25, RX23 = 100	C245445
RX26, RX24 = 100	C245445
CX1, CX4 = 10nF	C77053
CX2, CX5 = 10nF	C77053
CX3, CX6 = 10nF	C77053
RX21 = 1k	C2907372 C102060
RX4 = 1k	C2907372 C102060
RX11, RX15 = 120K	C163387
RX16, RX19 = 220K	C130275
RX28 = 180K	C137381
RX29 = 180K	C137381
D1 = 1N4148	C402212
D2 = 1N4148	C402212
CF = 1nF	C46553
U1, U2 = OPA192	C2861286
alternatively TLV07	C190218
CX = 680nF P=5mm Polyester Film Capacitor	C280298
Contact pins	C5199814
Male connector 8P (output) =	C225494

MYELINATED AXON





Figure 6. Schematics and PCB view (components side) of myelinated axon.

List of components of myelinated axon (with LCSC.com references):

R1 to R10 = 120K	C163387
R11 to R20 = 10M	C275643
RP1 = 1M15	C368017
RP2 = 120K	C163387
C1 to C10 = 220pF	C107081
Contact pins	C5199814
Female connector 5P (output)	C35167
Male connector 5P (input)	C138801

DEMYELINATED AXON





Figure 7. Schematics and PCB view (components side) of demyelinated axon.

List of components of myelinated axon (with LCSC.com references):

R1 to R10 = 120K

C163387

R11 to R20 = 1M

C108083

RP1 = not used

RP2 = not used

C1 to C10 = 10nF

C1589

Contact pins

C5199814

Female connector 5P (output)

C35167

Male connector 5P (input)

C138801

Important notes :

1. We suggest creating GERBER files immediately prior to PCB impression by importing the attached JSON files to EasyEDA.com. This allows also eventual modifications if needed.
2. Initially, the switches S1, TEA, PHR, TTX are in the LOW position and no predefined stimulation protocols are applied – the stimulation signal enters directly to the pin of the pipette on the SOMA PCB.
3. When the S1-switch is in HIGH position, the stimulation signal is deviated to the HDX4 connector. This allows stimulating the SOMA with another protocol via auxiliary PCB connected to HDX4. Attention: In the absence of this auxiliary PCB and with the switch S1 in HIGH position, the entry of the SOMA is at high impedance.
4. When S1 and PHR switches are in HIGH position, the pipette pin on SOMA PCB is disconnected. Instead, photoresistance (LDR, photoreceptor) becomes engaged. This allows observation of oscillatory behavior of SOMA under different illumination conditions.
5. S1 in LOW position and TTX in HIGH position, only passive responses of the SOMA are possible.
6. All switches in LOW position and with a subthreshold stimulus applied to the pipette – only passive response is observable.
7. PCB card of dendrites is connected to HDX2 socket. (for temporal and spatial summation experiments)
8. PCB cards of axones are connected to HDX3 socket. Node of Ranvier is connected in cascade to the axon.
9. The pin of output of the SOMA is referenced to AGND.